

Edit-R™ Gene Editing Workflows

Choose the right tools for your application

Whether your goal is a functional protein knockout for your gene of interest, running a loss-of-function screen, or creating an insertion or other knock-in, this workflow guide will assist you in selecting the right Edit-R™ gene editing reagents for your application.

Choose your application

Gene knockout

Gene knock-in

precise insertion or alteration of a gene

Single gene knockout

Knockout cell line creation or loss-of-function analysis in cell population

Choose a Cas9 reagent

Cas9 AAAA
Transfect Edit-R Cas9 mRNA or protein NLS for transient expression with no risk of DNA integration and fewer off-targets.

Cas9
Use Edit-R Lentiviral Cas9 nuclease for stable or inducible cell line creation for optimal editing efficiency.

Choose guide RNA

Synthetic guide RNAs work best for transient activity, and required no cloning or purification!


Optimize Cas9 delivery and expression AND guide RNA delivery

Cas9 + guide RNA transfected cells

3 days
Loss-of-function analysis in cell population.

3 days
Creating a knockout cell line? Single colony expansion in 96-well plates.

Assess gene editing efficiency and functional knockout phenotype



Loss-of-function screening of multiple genes at once

Pooled screening with no need for automation

Choose a Cas9 reagent and create stable cell line

Transduce lentiviral sgRNA pools

Perform pooled screen experiment

Analyze enriched sgRNA constructs in Reference vs. Experimental sample

Cas9
Choose lentiviral particles for the creation of stable or inducible Cas9 cell lines with choice of optimal promoter for your cell type.

sgRNA
Edit-R lentiviral sgRNA pools are fully sequenced libraries of algorithm-designed sgRNA demonstrating efficient gene editing at single-copy integrations.

Collect Reference sample, apply selective pressure to Experimental sample

Reference sample
Experimental sample
Experimental sample

Isolate gDNA
Reference gDNA
Experimental gDNA

Arrayed screening for one-gene-per-well analysis

Choose a Cas9 reagent and create stable cell line

Choose guide RNA and deliver to cells

Perform arrayed screen experiment

Assess loss-of-function phenotype

Cas9
Use Edit-R Lentiviral Cas9 nuclease for stable or inducible cell line creation for optimal editing efficiency.

Cas9
Edit-R synthetic guide RNA is available in arrayed predefined gene collections or custom libraries

Delivery with appropriate DharmaFECT transfection reagent formulation

3 days
Assess loss-of-function with a phenotypic assay

Choose a Cas9 reagent

Cas9 AAAA
Transfect Edit-R Cas9 mRNA or protein NLS for transient expression with no risk of DNA integration and fewer off-targets.

Choose guide RNA and donor oligo source

Synthetic guide RNAs work best for transient activity, and required no cloning or purification!

DNA donor plasmid for insertion of >50 nt
or
Single-stranded DNA oligo donor for insertions of <50 nt

Optimize Cas9 delivery and expression AND guide RNA as well as donor oligo delivery

Cas9 + custom guide RNA + donor DNA mixed population

Enrichment using FACS or antibiotic resistance (when applicable)
Expansion of selected cells

Assess HDR efficiency

3 days
Detection of HDR modification/insertion RFLP or junctional PCR

Single colony expansion into 96-well plates

2-3 weeks

Characterize clonal cell line

Ready to learn more about these offerings?

Contact your local Horizon Discovery representative or email us at ts.dharmacon@horizondiscovery.com